



DENTAL CERAMIC OVENS

zubler[®]
dental & technik

VARIO PRESS[®] 300.eZR



VARIO 200ZR



**FOR PERFECT PRESSINGS OF
LITHIUM DISILICATE CERAMICS**



**TTC MODE FOR
ACTIVE LINEAR COOLING**



Zubler Gerätebau GmbH
Buchbrunnenweg 26
D- 89081 Ulm

Tel.: + 49 (0) 731 - 14 52 0
Fax: + 49 (0) 731 - 14 52 13
www.zubler-group.de



zubler
group



Product comparison:



	V 200	V 200zR	VP 300	VP 300.e	VP 300.ezR
Ceramic oven	●	●	●	●	●
Press ceramic oven	○	○	●	●	●
ADVANCED PRESS™ process for lithium disilicate ceramics	○	○	○	●	●
Intelligent Press (i-press)	○	○	●	●	●
Z- Drying (pre-drying at constant temperature)	●	●	●	●	●
TTC - program for active, linear cooling (5°C/min - 45°C/min)	○	●	○	○	●
2-step programs with two free programmable heating rates	●	●	●	●	●
2-step programs with TTC cooling	○	●	○	○	●
PFC (Power Fail Control)	●	●	●	●	●
500 programs, each individually programmable	●	●	●	●	●
USB connection	●	●	●	●	●
Impact and heat-resistant LCD glass touch display	●	●	●	●	●
Custom pre-programmed for the client	○	●	○	●	●
Accessory set (honeycomb carrier, pins)	○	●	○	●	●
2-year warranty	●	●	●	●	●

