

Make good press results perfect press results!

**The intelligent power pack is the problem solver
for malpositions of ceramic objects**

- Immediately after opening the press furnace place the pressed muffle in the **VARIO SC Surface Conditioner**.
- Compensation of solidification contractions and thereby avoiding misalignments in the pressed object.
- Defined, even cooling of the pressed muffle at reduction of the cooling time by approx. 50%



Best order also now...

CONCEPT VEST™ - the special investment for lithium disilicate ceramics

The investment especially developed for processing of lithium disilicate glass ceramics guarantees ...

- perfect surfaces
- a minimal reaction layer after pressing
- detailed reproduction of the pressed objects
- consistent expansion values through strictest quality controls



25 x 100 g

Item-No. 898/2066

50 x 100 g

Item-No. 898/2065

Mixing Liquid / 1 Liter

Item-No. 898/204

Zubler® Disposable Plungers

For pressing of all common dental ceramics and **lithium disilicate**.

Their excellent quality and the exact dimensions guarantee the repeatable success with pressings in muffles of all sizes.



ø 12 mm

Item-No. 20/0235

ø 13 mm

Item-No. 20/0238

ø 16 mm

Item-No. 20/0237

Zubler® Investment System for Lithium Disilicate Ceramics

The ribbed surface of the muffles ensures faster heat absorption and release. Sector lines make it easier to locate the pressed object, damage during devesting is almost impossible. The Flexing is easy to clean and has a long service life.



Flexring 100 g Set

Item-No. 20/0228

Flexring 200 g Set

Item-No. 20/0229

